

PRESTAN Professional Jaw Thrust Manikin

Frequently Asked Questions



What is different about the Jaw Thrust Manikin?

Trainers who teach the jaw thrust method will benefit from this patented feature. The uniqueness of the PRESTAN Professional Jaw Thrust Manikin is in the design of the jawline, allowing the jaw to be thrust forward to open the airway. This ability is in addition to the head tilt. Available in Adult size only.

Can I purchase the Jaw Thrust Head for my Manikin?

Yes. In addition to purchasing the PRESTAN Professional Jaw Thrust Manikin, the Jaw Thrust Heads may be purchased separately. They are interchangeable with the standard PRESTAN manikin heads for the PRESTAN Professional Adult Manikins.

Do the Jaw Thrust Manikins comply with the industry guidelines?

The Jaw Thrust Heads are on Professional Adult Manikin torsos. Therefore, they are compliant with the industry guidelines if they have the CPR Rate Monitor installed.